SUPER ERECTA SHELF®

Inhalation Therapy Cart

□ Super Erecta Shelf® Construction: Allows custom height adjustment and positive locking without the use of nuts, bolts or setscrews and without the use of special tools.

Note: To maintain the conductivity of the cart, aluminum split sleeves must be used rather than plastic.

□ PVC-Coated Top & Bottom Inserts: The compartmented inserts, which snap over the top and bottom frames, have a resilient PVC coating which aids in reducing vibration and noise.

☐ Chrome-Plated Finish on handles, frames and shelf.


□ Accommodates Any Combination of twelve D & E cylinders. Each PVC-coated insert has twelve compartments, each measuring 41/2" (114mm) square.

5" (127mm) Conductive Casters. See caster NOTICE on assembly instruction sheet.

Copyright © 1984 InterMetro Industries Corp.



ITC12C



DIMENSIONS

	Width	Length	Height	Approx. Pkd, Wt.
Cat. No.	(in.) (mm)	(in.) (mm)	(in.) (mm)	(lbs.) (kg)
ITC12C	201/2 520	26 660	39 990	401/2 18